


WHY WAIT FOR 4K, IT'S HERE!

NATIONWIDE COVERAGE FOR PRODUCT INSTALLATION, WARRANTY, AND SERVICE




*Please contact Prescott's for information
regarding your local representative.*

"The Prescott's OmniPlus/OSM200 is an excellent office microscope for the ENT surgeon. It is very affordably priced and functions exceptionally well. I give Prescott's my highest recommendation." ~Dr. Kuperan, Houston Physicians Hospital

"Good clear optics at a good price!" ~Dr. Andrew Diamond, WMG Northside ENT


"Scope function is easy and simple. Visual acuity is top notch, and LED lighting is exactly what we needed." ~ Dr. Montano, McAllen Surgical Specialty Center



OmniPlus/OSM200



Prescott's Microscopes
1770 Deer Creek Rd.
Monument, CO 80132
844-891-1897



EMERGO EUROPE
Prinsesdag 20
2514 AP The Hague
The Netherlands

PRESCOTT'S, INC.
Surgical Microscope Specialists


PRESCOTT'S, INC.
Surgical Microscope Specialists

OmniPlus/OSM200

Incorporating everything needed to fully concentrate on the procedure at hand. For better results. Everyday.

With its compact, flexible and robust design, the OSM200 stand is suited for any size operating room.

Magnification	5 Step Apochromatic magnification changer (0.4, 0.6, 1.0, 1.6, 2.5)
Focus	Fine focus objective lens
Main Lens	f= 250mm, Apochromatic
Illumination	60 watt LED (5500 K) 12v 150 Halogen bulb
Binocular	f= 170mm, 10x straight binocular
Balancing	Mechanical
Asepsis	Drapes, sterilizable caps
Suspension	Height: 55 inches Base: 23.5" wide, 24" deep
Video	4K Recording and still image capture
Accessories	"Endo Wedge" for use in Endodontics Video Remote, 1-180 Inclinable Binocular, Variable objective lens: 207mm-407mm



LED Illumination

Superior brightness approximating natural light, CRI 92
Long lasting bulb, estimated life 50,000 hours
5700K Color Temperature
Turret Design on LED will accept multiple cable sizes
Silicon based liquid fiber optic cable for improved light




Focus/Magnification

Fine focus objective lens
12 mm of focus travel
Multiple working distances available
5 Step Apochromatic mag changer
Real Magnification 2.7X through 17X




Endodontics

Combining the Prescott's "Endo Wedge" with inclinable 0-180 binocular provides surgeon with maximum flexibility in dental



Wall Mount

When procedure space is at a premium take advantage of Prescott's wall mount, while sacrificing nothing in performance.



Video

Integrated 4K Streaming video and still image capturing
64G memory
Integrated camera can be paired with phone, tablet, or remote control
Designated HDMI cable for monitor display
USB 3.0/SD Card Reader